

LEYSPEC Residual Gas Analyzer (RGA)

Quick Start Manual 300781309_002_C0

Part Numbers

220010V01
220011V01
220012V01
220013V01
220020V01
220021V01



Contents

	Page
0 Introduction	5
0.1 Before use	5
1 Equipment configuration	6
1.1 Standard equipment configuration	6
2 Part name and component function	7
2.1 Operation switches and LEDs	7
2.2 Rear panel	8
2.3 DIP switch	8
3 Installation and connection	11
3.1 Preliminary operation	11
3.2 Installation	11
4 Measurement	16
4.1 Using a computer	16
4.2 Using the sensor unit	17
5 Warranty	18
5.1 Items covered under warranty	18
5.2 Duration of guarantee	18
5.3 Scope	18
5.4 Failures which void the warranty	18
5.5 Other	18
6 Certificate of decontamination	19

Introduction

0 Introduction

This manual is to be used as a quick check of the operation and display of the Residual Gas Analyzer.

Refer to the instruction manual located on the CD for detailed information on proper operation, precautions and safe use of the unit.

0.1 Before use

Confirm that the PC has a LAN port. To use a USB/LAN converter, install a driver.

The sensor unit is a 24 V d.c. drive. Prepare the AC/DC adaptor or 24 V a.c./d.c. power supply (more than 50 W).

Confirm the model to be used. The model can be found on the sensor connector panel side of the sensor unit and flange.

Equipment configuration

1 Equipment configuration


1.1 Standard equipment configuration

Item	Quantity
Sensor unit	1 piece
Sensor tube	1 piece
Software and instruction manual CD	1 piece
AC/DC adaptor	1 piece
LAN cable	1 piece

Part name and component function

2 Part name and component function


2.1 Operation switches and LEDs



Item	Description	
1	Analyzer tube connector	For connection to the analyzer tube.
2	Measurement value/fault display	Shows measured partial pressure value. For faults, displays the error code.
3	Gas species/total pressure display LED	Displays the gas species and total pressure being measured.
4	POWER	POWER LED (green) comes on when power is applied.
	SEM	Secondary Electron Multiplier (SEM) LED (green) comes on when SEM is used.
	ERROR	ERROR LED (red) comes on if a fault occurs.
5	Gas species/total pressure selection button (CH)	Select the gas species you want to measure or total pressure. The orange LED shows the current selection.
6	Measurement start button (START)	Starts measurement. The upper orange LED indicates that measurement has started.
7	Measurement value units display LED	The unit of measured partial pressure. For O2/N2 the display shows a percentage (%).

Part name and component function


2.2 Rear panel



Item	Name	Function
1	Power fuse holder (FUSE)	Over-current protection (Time lag type 6.3 A).
2	Rotary switch for regulating tuning voltage (TUNING)	Regulates RF voltage tuning with the analyzer tube.
3	Tuning voltage check terminal (+, -)	Terminal for checking RF voltage tuning with the analyzer tube.
4	Power inlet connector (POWER IN)	Power cable connector (includes ground).
5	Resolution adjustment potentiometer (RESOLUTION L, H)	L: Adjusts the resolution of the low mass number. H: Adjusts the resolution of the high mass number.
6	External input/output connector (EXT-I/O)	Signal I/O connector (D-sub 15-pin male).
7	Initial setting DIP switch (CONFIG)	DIP switch for SEM selection, IP address setting, the unit address setting and local/remote selection.
8	Ethernet connector (LAN)	RJ45 connector (8-pole 8-wire) for connecting the LAN cable

2.3 DIP switch

Set up the DIP switch numbers 1, 3, and 4 through 8 with the power supply OFF.



Part name and component function

No.	Function
1	To change IP address between "Fixed IP address" or "Automatic IP address". ON = automatic OFF = fixed
2	Can be changed with power on and measurement stopped. To select between secondary electron multiplier (SEM) or Faraday Cup (FC). ON = SEM OFF = FC
3	Not used.
4	Can be changed with power on and measurement stopped. To set local and remote operation. ON = Local mode (control panel operation only). OFF = Remote mode (PC operation)
5	
6	To set the sensor address for this LEYSPEC Residual Gas Analysers unit.
7	See the table below for switch settings.
8	

Address	S5	S6	S7	S8	Address	S5	S6	S7	S8
1	0	0	0	1	9	1	0	0	1
2	0	0	1	0	10	1	0	1	0
3	0	0	1	1	11	1	0	1	1
4	0	1	0	0	12	1	1	0	0
5	0	1	0	1	13	1	1	0	1
6	0	1	1	0	14	1	1	1	0
7	0	1	1	1	15	1	1	1	1
8	1	0	0	0	16	0	0	0	0

The following table shows the settings before shipment from the factory.

LEYSPEC view100S/200S & LEYSPEC ultra 200S/300S

S1	S2	S3	S4	S5	S6	S7	S8
OFF	ON	OFF	OFF	OFF	OFF	OFF	ON

LEYSPEC view100/200

S1	S2	S3	S4	S5	S6	S7	S8
OFF	OFF	OFF	OFF	OFF	OFF	OFF	ON

Part name and component function


IP address setting	: OFF (fixed)
Local/Remote	: REMOTE
Address	: 1

List of fixed IP addresses:

Sensor1	: 192.168.250.11	Sensor9	: 192.168.250.19
Sensor2	: 192.168.250.12	Sensor10	: 192.168.250.20
Sensor3	: 192.168.250.13	Sensor11	: 192.168.250.21
Sensor4	: 192.168.250.14	Sensor12	: 192.168.250.22
Sensor5	: 192.168.250.15	Sensor13	: 192.168.250.23
Sensor6	: 192.168.250.16	Sensor14	: 192.168.250.24
Sensor7	: 192.168.250.17	Sensor15	: 192.168.250.25
Sensor8	: 192.168.250.18	Sensor16	: 192.168.250.26

Installation and connection

3 Installation and connection



* This figure is for the single sensor type.

- A. LEYSPEC Residual Gas Analysers
- B. LAN cable
- C. AC/DC adaptor
- D. PC
- E. Mains power

3.1 Preliminary operation

1. Unpack the instrument and check quantities.
2. Check if any components are damaged.

3.2 Installation

3.2.1 Install the sensor

1. Measuring position
 - This product measures pressure by measuring the static pressure in the position where the analyzer tube is connected. Take care in selecting the measuring position because the measurement value is affected if there is a flow in the vacuum system, source of out-gas or source of high intensity electrons or ions.
 - The correct pressure measurement cannot be made or the analyzer tube may fail if the analyzer tube is subject to vibration, heat radiation, high intensity electromagnetic field or high intensity radiation.
 - The LEYSPEC Residual Gas Analysers must be connected to power supplies or instruments that conform to the requirements of a grounded protective extra-low voltage (SELV-E according to EN 61010).
 - Keep enough room around the unit as to not cover panel vents and allow for sufficient heat dissipation.
2. Installing the analyzer tube

Installation and connection

- Install the analyzer tube so that the aperture plane of the analyzer tube is installed in parallel to the flow of gas. Especially, see to it that gas is not introduced into the analyzer tube in a beam form. Also make sure that foreign objects cannot fall into the analyzer tube by installing a mesh if the aperture faces up.
- If you touch the gasket with bare hands or use a contaminated gasket it can cause measurement errors. Use a new gasket and, before use, wipe it with a clean cloth wetted with alcohol. To avoid leaks, never reuse gaskets.
- If the filament is subjected to strong lateral impact or vibration when it is hot, it may break or contact the grid electrode.
- If the analyzer tube is connected at the end of a long or thin pipe, the pressure in the analyzer tube will increase because of the reduced conductance. This can cause measurement errors or filament burnout.
- As indicated in [Operation switches and LEDs](#) on page 3, fix the analytical tube after paying attention to the orientation of the operation panel on the sensor unit. To assist you, a "Display Side" label is attached to the flange of the analyzer tube.

CAUTION:

After fixing the analysis tube, do not attempt to twist it in order to change the orientation of the sensor unit. This can cause the sensor unit to malfunction.



3.2.2 Install the sensor unit

CAUTION:

Ensure the sensor unit does not strike other structures.



CAUTION:

Attach or detach the sensor unit from the analyzer tube using the knurled ring. Do not rotate the sensor unit, this can break the connector or the wiring.




CAUTION:

Only ever tighten the knurled ring by hand. If you use tools, there is a risk of making it impossible to detach, breaking the connector, or damaging the internal wiring.



Installation and connection

1. Attach the sensor unit and the analyzer tube and tighten the knurled ring by hand.



3.2.3 Install the power connector




WARNING:

Make connections to the correct pins when wiring the power connector. If you apply power to the incorrect pins internal damage could result.

1. Wire the power connector as follows:
 - Power +24 V
Pin that supplies 24 V d.c. power.
 - Power GND
GND when +24 V d.c. power is supplied.
 - Frame GND
Frame ground. Connect it to the instrument frame.
Connect the earth connection to Class A ground to prevent leakage currents.


Installation and connection



Item	Function	Description
1	Frame GND	Frame ground. Connect it to the instrument frame. (AC adaptor code; Black + Green dot)
2	Power GND	GND when +24 Vd.c. power is supplied. (AC adaptor code; black)
3	Power +24 V	Pin that supplies 24 Vd.c. power. (AC adaptor code; white)

3.2.4 LAN cable connection

To connect a single sensor, select LAN. Connect the cable, as shown in the figure below.




1. PC
2. LAN connector
3. LAN cable
4. Sensor unit

Installation and connection

3.2.5 LAN cable connection

To connect multiple sensors, select LAN. Set the address and LAN terminator. Connect the cable.



1. PC
2. LAN connector
3. LAN cable
4. Sensor unit
5. Switching hub
6. Sensor unit (maximum 16 units)

CS/1871/A

4 Measurement

There are 2 measurement methods you can use on the Residual Gas Analyser:

1. Connect the unit to a computer (PC) in which the LEYSPEC software gas analysis software is installed.
2. Use the push button, indicator light and partial pressure display value on the sensor unit.

4.1 Using a computer

1. Set the PC's IP address. When using the fixed IP address of the sensor unit (192.168.250.11 in case of Sensor 1), open the property of the IP address (TCP / IP-Ver4) from the network setting of the PC and set the IP address = 192.168.250.xxx (xxx is a numerical value other than "11" of the sensor unit). Refer to the instruction manual for details such as setting variable IP address.

- When connecting to a PC, confirm the setting of both IP addresses.
- Turn off the Windows Firewall.
- Communication can only be used with a wired connection with Windows PC and the dedicated software installed. Communication with communication devices other than a Windows PC is not supported.

2. Install the LEYSPEC software (called LEYSPEC Software in this manual) gas analysis software on your PC.

3. Use a LAN cable to connect the LAN port (RJ45) of the PC to the LAN connector of the sensor unit.

4. After supplying power to the LEYSPEC Residual Gas Analysers unit, double-click the LEYSPEC Software icon on the desktop of the PC and run the LEYSPEC Software software.

5. Enter the sensor number and type in the sensor type dialog, displayed first, and press the "OK" button.

6. When communication between the PC and the LEYSPEC Residual Gas Analysers is established, LEYSPEC Software loads the default recipe or the recipe used last.

7. Press the "FIL" button to turn on the filament of the analyzer tube. The "F" and "R" lights come on.

- Before pressing the START button, check the operating pressure (about 1×10^{-3} Pa) on a separate gauge.
- If the filament is turned on too soon after pumping, the high pressure protection of the LEYSPEC Residual Gas Analysers might work because of the emitted gas of the filament. In that case, confirm operating pressure again and turn on the filament several times to reduce the emitted gas.

Measurement

8. Press the "MES" button to start gas analysis. The "S" and "M" lights come on. (Whether the F, R and S lights are on or off depends on the settings on the control panel.)



9. To terminate gas analysis, press the "MES" button again.

10. Turn off the "FIL" button and then turn off the analyzer tube filament.

11. Exit the software when measurement is completed.

4.2 Using the sensor unit

Refer to the operation manual supplied on the software installation CD.

5 Warranty

This product was shipped after company inspection. However, if a failure occurs such as a defect in manufacturing or damage during transportation, contact your supplier.

5.1 Items covered under warranty

This unit.

5.2 Duration of guarantee

One (1) year after shipping date.

5.3 Scope

1. Damage to product caused by a failure on delivery.
2. Product not meeting the standard specifications when used under normal service conditions as stated in the manual.

5.4 Failures which void the warranty

1. Failure occurred after expiration of the warranty period.
2. Failure caused by force of nature, such as fire, storm and flood damage, earthquake, lightning strike, and so forth.
3. Failure occurred due to careless handling or incorrect usage.
4. Product is refurbished, disassembled or repaired without supplier acceptance.
5. Failure occurred in an abnormal environment, such as an intense electromagnetic field, radiation, high temperature, high humidity, flammable gases, corrosive gases, dust, and so forth.
6. Failure due to noise.
7. Ion source and secondary electron multiplier and sensor being used (expiration of life, measurement error, and so forth).
8. Cable in use (cable burnout due to improper installation, poor contact, and so forth).

5.5 Other

1. For questions and consultation, the customer can check the model and serial number and ask the local supplier representative.
2. The content of this document is subject to change without notice.

Certificate of decontamination

6 Certificate of decontamination

All material must be certified as decontaminated and this certificate must be submitted to the closest supplier service centre or sales office before shipment. Refer to the CD.



WARNING:

Read this instruction manual before installing, operating, inspecting, or servicing this product to familiarize yourself with safety precautions, specifications and operations.

Sales and Service

Germany

Leybold GmbH

Bonner Strasse 498
D-50968 Cologne
T: +49-(0)221-347 1234
F: +49-(0)221-347 31234
sales@leybold.com
www.leybold.com

Leybold GmbH

Sales Area North
Branch Office Berlin
Industriestrasse 10b
D-12099 Berlin
T: +49-(0)30-435 609 0
F: +49-(0)30-435 609 10
sales.bn@leybold.com

Leybold GmbH

Sales Office South
Branch Office Munich
Karl-Hammerschmidt-Strasse 34
D-85609 Aschheim-Dornach
T: +49-(0)89-357 33 9-10
F: +49-(0)89-357 33 9-33
sales.mn@leybold.com
service.mn@leybold.com

Leybold Dresden GmbH

Service Competence Center
Zur Wetterwarte 50, Haus 304
D-01109 Dresden
Service:
T: +49-(0)351-88 55 00
F: +49-(0)351-88 55 041
info.dr@leybold.com

Europe

Belgium

Leybold Nederland B.V.

Belgisch bijkantoor
Leuvensesteenweg 542-9A
B-1930 Zaventem
Sales:
T: +32-2-711 00 83
F: +32-2-720 83 38
sales.zv@leybold.com
Service:
T: +32-2-711 00 82
F: +32-2-720 83 38
service.zv@leybold.com

France

Leybold France S.A.S.

Parc du Technopolis, Bâtiment Beta
3, Avenue du Canada
F-91940 Les Ulis cedex
Sales and Service:
T: +33-1-69 82 48 00
F: +33-1-69 07 57 38
info.ctb@leybold.com
sales.ctb@leybold.com

Leybold France S.A.S.

Valence Factory
640, Rue A. Bergès
B.P. 107
F-26501 Bourg-lès-Valence Cedex
T: +33-4-75 82 33 00
F: +33-4-75 82 92 69
marketing.vc@leybold.com

Great Britain

Leybold UK LTD.

Unit 9
Silverglade Business Park
Leatherhead Road
Chessington
Surrey (London)
KT9 2QL
Sales:
T: +44-13-7273 7300
F: +44-13-7273 7301
sales.ln@leybold.com
Service:
T: +44-13-7273 7320
F: +44-13-7273 7303
service.ln@leybold.com

Italy

Leybold Italia S.r.l.

Via Filippo Brunelleschi 2
I-20093 Cologno Monzese
Sales:
T: +39-02-27 22 31
F: +39-02-27 20 96 41
sales.mi@leybold.com
Service:
T: +39-02-27 22 31
F: +39-02-27 22 32 17
service.mi@leybold.com

Netherlands

Leybold Nederland B.V.

Floridadreef 102
NL-3565 AM Utrecht
Sales and Service:
T: +31-(30) 242 63 30
F: +31-(30) 242 63 31
sales.ut@leybold.com
service.ut@leybold.com

Switzerland

Leybold Schweiz AG

Hinterbergstrasse 56
CH-6312 Steinhausen
Warehouse and shipping address:
Riedthofstrasse 214
CH-8105 Regensdorf
Sales:
T: +41-44-308 40 50
F: +41-44-302 43 73
sales.zh@leybold.com
Service:
T: +41-44-308 40 62
F: +41-44-308 40 60
service.zh@leybold.com

Spain

Leybold Hispánica, S.A.

C/. Huelva, 7
E-08940 Cornellà de Llobregat
(Barcelona)
Sales:
T: +34-93-666 43 11
F: +34-93-666 43 70
sales.ba@leybold.com
Service:
T: +34-93-666 46 11
F: +34-93-685 43 70
service.ba@leybold.com

America

USA

Leybold USA Inc.

5700 Mellon Road
USA-Export, PA 15632
T: +1-724-327-5700
F: +1-724-325-3577
info.ex@leybold.com
Sales:
T: +1-724-327-5700
F: +1-724-333-1217
Service:
T: +1-724-327-5700
F: +1-724-325-3577

Brazil

Leybold do Brasil

Rod. Vice-Prefeito Hermenegildo Tonolli,
nº. 4413 - 6B
Distrito Industrial
Jundiaí - SP
CEP 13.213-086
Sales and Service:
T: +55 11 3395 3180
F: +55 11 99467 5934
sales.ju@leybold.com
service.ju@leybold.com

Asia

P. R. China

Leybold (Tianjin)

International Trade Co. Ltd.

Beichen Economic
Development Area (BEDA),
No. 8 Western Shuangchen Road
Tianjin 300400
China
Sales and Service:
T: +86-400 038 8989
T: +86-800 818 0033
F: +86-22-2697 4061
F: +86-22-2697 2017
sales.tj@leybold.com
service.tj@leybold.com

India

Leybold India Pvt Ltd.

T-97/2, MIDC Bhosari
Pune-411 026
Indien
Sales and Service:
T: +91-80-2783 9925
F: +91-80-2783 9926
sales.bgl@leybold.com
service.bgl@leybold.com

Japan

Leybold Japan Co., Ltd.

Headquarters
Shin-Yokohama A.K.Bldg., 4th floor
3-23-3, Shin-Yokohama
Kohoku-ku, Yokohama-shi
Kanagawa-ken 222-0033
Japan
Sales:
T: +81-45-471-3330
F: +81-45-471-3323
sales.yh@leybold.com

Leybold Japan Co., Ltd.

Tsukuba Technical Service Center
1959, Kami-yokoba
Tsukuba-shi, Ibaraki-shi 305-0854
Japan
Service:
T: +81-29 839 5480
F: +81-29 839 5485
service.iik@leybold.com

Malaysia

Leybold Malaysia

Leybold Singapore Pte Ltd.

No. 1 Jalan Hi-Tech 2/6
Kulim Hi-Tech Park
Kulim, Kedah Darul
Aman 09000
Malaysia
Sales and Service:
T: +604 4020 222
F: +604 4020 221
sales.ku@leybold.com
service.ku@leybold.com

South Korea

Leybold Korea Ltd.

3F, Jellzone 2 Tower
Jeongja-dong 159-4
Bundang-gu Sungnam-si
Gyeonggi-do
Bundang 463-384, Korea
Sales:
T: +82-31 785 1367
F: +82-31 785 1359
sales.bd@leybold.com
Service:
623-7, Upsung-Dong
Cheonan-Si
Chungcheongnam-Do
Korea 330-290
T: +82-41 589 3035
F: +82-41 588 0166
service.cn@leybold.com

Singapore

Leybold Singapore Pte Ltd.

42 Loyang Drive
Loyang Industrial Estate
Singapore 508962
Singapore
Sales and Service:
T: +65-6303 7030
F: +65-6773 0039
sales.sg@leybold.com
service.sg@leybold.com

Taiwan

Leybold Taiwan Ltd.

10F., No. 32, Chenggong 12th St.,
Zhubei City, Hsinchu County 302
Taiwan, R.O.C.
Sales and Service:
T: +886-3-500 1688
F: +886-3-550 6523
sales.hc@leybold.com
service.hc@leybold.com

Headquarter

Leybold GmbH

Bonner Strasse 498
D-50968 Cologne
T: +49-(0)221-347-0
F: +49-(0)221-347-1250
info@leybold.com



www.leybold.com